



virus safety in plasma- derived therapeutics:

A MERCK PERSPECTIVE

IPFA YOGYAKARTA, March 2017

Louis Wong

Associate director, Plasma Initiative, Asia-Pacific



MERCK

Content


- **Regulatory Expectations for Plasma-Derived Therapeutics**
 - **Multi-Layer Strategy for Virus Safety**
- Prevention - Raw Materials and Chemicals**
- Removal – Clearance Technologies**
- Detect – In-Process and Final Product Testing**
- **BioReliance Testing Services**



Virus Contamination in Plasma – A Reality

Plasma Product Contamination

- 1980s: HIV contaminated Coagulation Factor
- 1994: Hepatitis C contaminated IgG
- 1995: vCJD contamination
- Other threats more recently... Hepatitis E, West Nile Virus, Ebola, Zika...



Virus Safety in Plasma : Part of a Total Risk Mitigation Program

- ✓ Regular, qualified donors
- ✓ Donation control (serological and NAT tests)
- ✓ Inventory Hold, look back procedure
- ✓ Serological and NAT Tests on plasma pool
- ✓ **cGMP**
- ✓ **Pathogen inactivation/removal steps**
- ✓ **Batch Release**
- ✓ Post-marketing surveillance, pharmacovigilance





regulatory expectations for human plasma protein products



MERCK

Regulatory Expectations for Human Plasma Protein Products

- Guideline on plasma-derived medicinal products. EMA/CHMP/BWP/706271/2010 (2011)
- ✓ Incorporate at least 2 orthogonal, effective virus inactivation/removal steps with at least 4 log virus reduction in each step is desirable
- ✓ Virus Inactivation methods/Virus removing filters

Virus Group Choices	Types	Details
Enveloped Viruses	HIV 1 Model for Hepatitis C	Pestivirus BVDV could be considered a “worst-case” for HCV.
Enveloped DNA Viruses	Herpes Virus Hepatitis B (duck model)	A validation study should be performed with an appropriate enveloped DNA virus, e.g. a herpesvirus such as pseudorabies.
Non-enveloped Viruses	Hepatitis A B19V	HAV and B19V should be used for validation studies for coagulation factors (History of transmission).
Model viruses for virus reduction filtration (nanofiltration)	HIV BVDV	For small pore size filters designed for removal of small non-enveloped viruses, HIV and BVDV should still be part of the virus panel, but robustness studies may focus on small non enveloped viruses. For medium pore size filters, BVDV is appropriate for robustness studies.



Virus inactivation/removal capacity

Points to consider for specific product classes

EMA/CHMP/BWP/706271/2010 (2011) / Section 8.3

Product Class	Details
Coagulation Factors	Non-enveloped viruses such as hepatitis A and B19V have been transmitted by this class of products. For Factor IX products, steps should be included in the process that are effective for HAV and B19V.
Immunoglobulins	Immunoglobulin products have a good safety record for the known non-enveloped viruses due in part to the contribution from neutralising antibodies in the product.
Albumin	Albumin manufactured by an established fractionation process that includes the terminal pasteurisation specified in the European Pharmacopoeia monograph, has an excellent virus safety record.
S/D Plasma	SD plasma has good safety measures for enveloped viruses and safety measures are in place for HAV and B19 (Ph. Eur. monograph Human Plasma (Pooled and Treated for Virus Inactivation)).



Manufacturing Process Validation Activities by Drug Development Phase

Process Validation Activities	Clinical Phases	Supporting Activities
Perform Preliminary Clearance of Process Impurities	1	---- Initiate Process and Assay Development Initiate Preliminary Process Characterization Establish Preliminary Reference Standard
---- Initiate Virus Clearance Studies		---- Initiate Clinical Manufacturing
----		---- Initiate Stability Studies
		Submit IND Application
Initiate Manufacturing Process Characterization Studies (test limits)	2	Conduct Continuing Process & Assay Development
Approve Manufacturing Process Validation Plan		Qualify or Validate Assays
Develop/Approve Process Validation Protocols		---- Establish Gene Therapy Product Characterization ---- Establish Manufacturing Process Definition
Conduct Additional Virus Clearance Studies		---- Conduct Pivotal Clinical Trails



Virus Safety: Multi-Layer Strategy A Merck Perspective

1



Prevent

**Virus Safety of Raw
Materials**

2



Detect

**Testing of In-Process
and Final Product**

3



Remove


**Virus Clearance
Technologies**



Overcoming challenges to transition your raw materials

Start with the end in mind

- *Quality standards differ for raw materials*
- *Quality requirements increase for drugs approaching commercial launch*





Overcoming challenges to transition your raw materials

Comparison of approaches

Traditional development

- ✓ Focus on end product safety and performance
- ✓ Fixed process – do not change make any changes after Research
- ✓ Fixed specifications
- ✓ Work with suppliers starting phase II



Enhanced development

- ✓ Systematic – understand material features and process to meet critical to quality attributes
- ✓ Supply Chain Risk Management to ensure uninterrupted drug supply
- ✓ Work with suppliers who can support you in transitioning to Clinical and comply with Regulatory requirements



The Merck EMPROVE® Concept


The EMPROVE concept brings a package of benefits to the plasma fractionation customer that include:

- Strong process safety and reliability – High quality raw materials with GMP production;
- Strong regulatory support – EMPROVE dossier in CTD format allows the customer to save time and money;
- Full and permanent access to main documents in our website;
- Top level Quality Control procedures, covering the manufacturing and the supply chain steps.

Quality, Safety and Stability are the key words in the EMPROVE® brand products for blood plasma fractionation



General Fractionation Process Steps with EMPROVE® Chemicals



Virus Safety: A Multi-Layer Strategy



Prevent

Virus Safety of Raw Materials



Detect

Testing of In-Process and Final Product



Remove

Virus Clearance Technologies



Types of Process Steps

EMA/CHMP/BWP/706271/2010 Guideline on plasma-derived medicinal products

“For all plasma-derived medicinal products, it is an objective to incorporate effective steps for inactivation/removal of a wide range of viruses of different physico-chemical characteristics.”

Inactivation

- Heat Treatment
 - Pasteurization (60°C, 10 hours, Albumin)
 - Lyophilization / dry heat
- Solvent/Detergent (TnBP, Triton X-100, Tween 80)
- Low/High pH (Caprylic acid, Acetic acid)
 - Column Elution/Sanitization

Partitioning

- Precipitation
 - EtOH
 - PEG
 - Caprylic Acid
- Chromatography
 - Ion Exchange
 - Affinity
 - Gel Filtration
- NanoFiltration



Viral Reduction Capacities by Contributive Steps

Spiked viruses		ENVELOPED VIRUSES				NON-ENVELOPED VIRUSES			
		HIV-1	Sindbis	BVDV	PRV	SV-40	EMCV	HAV	PPV
Model for		HIV	HCV		Enveloped DNA virus (HBV)	Highly resistant virus	HAV		Parvovirus B19
Specific steps	S/D treatment	e 4.4	e 5.4	NT	e 4.3	N/A	N/A	N/A	N/A
	Nanofiltration through 20 nm filter	5.6 [Ⓞ]	5.6 [Ⓞ]	5.6	5.6 [Ⓞ]	5.4 [Ⓞ]	e 5.2	NT	4.3
Contributive steps	Caprylic acid fractionation	e 4.0	NT	5.1	e 5.0	NT	e 5.6	e 5.6	3.7
	Anion-exchange chromatography	NT	NT	NT	NT	NT	1.3	NT	3.8
	Final product: Low pH incubation in the final container	4.0	NT	3.2	NT	NT	NT	NT	NT
	Final product: Virus-neutralising capacity by the antibodies	N/A	N/A	N/A	N/A	N/A	N/A	3.3 [Ⓞ]	N/A
Overall viral reduction capacity		e 18.0	e 11.0	13.9	e 14.9	5.4	e 12.1	e 8.9	11.8


Source: Production of ClairYg, Christophe Segard, LFB



Solvent/Detergent Viral Inactivation using Single-Use Technology



Objectives:

- Quantify Impact of S/D Treatment **on Chemical Compatibility, Extractables/Leachables** and **Mixing Efficiency** in Single-Use Bags
- Determine **Efficacy** of S/D Viral Inactivation in Single-Use Bags
- Investigate Impact on **Protein Activity**
- Create **Best Practices** for performing S/D Process Step in Single-Use Mixers




Outcomes of Study

Virus Reduction on Human IgG performed in Mobius SU bags



Mixing Efficiency in Mobius SU-Mixer 50L and 500L



Virus Removal by Nanofiltration

- Virus Filtration performed on 10 different Grifols plasma proteins/products
- Conclusion of Nanofiltration being extremely effective step in virus reduction
- Minimal influence of process parameters like pH, temperature, Conductivity, during NF

Protein (or product)	Enveloped viruses (≥ 40 nm)				Non-enveloped viruses		
					Small (< 35 nm)	Very small (< 25 nm)	
	PRV	HIV-1	WNV	BVDV	HAV/EMCV	PPV	
A	≥ 4.6	≥ 4.8	≥ 3.6	≥ 4.7	≥ 5.9	4.6	
B	≥ 6.1	≥ 5.6	≥ 4.5	≥ 4.5	5.2	4.4	
C	≥ 6.0	≥ 4.0	≥ 5.4	≥ 4.9	6.6	6.1	
D	D1	≥ 5.5	≥ 5.9	≥ 5.8	≥ 4.1	4.1	
	D2	≥ 5.2	≥ 6.9	n.d.	≥ 4.7	≥ 5.5	4.2
E	E1	≥ 4.2	≥ 3.8	≥ 6.2	≥ 3.7	≥ 6.4	5.0
	E2	≥ 5.4	≥ 5.4	n.d.	≥ 4.8	≥ 4.2	6.5
F	F1	≥ 6.0	≥ 6.8	≥ 7.0	≥ 6.3	≥ 5.3	3.9
	F2	n.d.	≥ 4.5	≥ 6.0	4.6	≥ 4.4	3.8
G	n.d.	≥ 6.3	n.d.	≥ 5.1	≥ 5.6	4.6	

BVDV: bovine viral diarrhea virus; EMCV: murine encephalomyocarditis virus; HAV: hepatitis A virus; HIV: human immunodeficiency virus; n.d.: not determined; PPV: porcine parvovirus; PRV: pseudorabies virus; WNV: West Nile virus.



Virus Clearance Validation



Prevent

Virus Safety of Raw Materials



Detect

Testing of In-Process and Final Product



Remove

Virus Clearance Technologies



BioReliance[®]


by **SAFC**



BioReliance Clearance Services

Perform studies at 2 separate facilities

- Dedicated facilities in both Rockville (MD), USA and Stirling, Scotland.
 - Rockville facility has 4 study suites and a BSL-3 lab.
 - Stirling facility has 5 study suites and a BSL-3 lab.
 - Both make use of state-of-the-art equipment Endpoint testing conducted in adjacent virus titration labs to increase the speed and efficiency of sample analysis.
- Trained and experienced personnel provide technical and regulatory support for all projects.
- Besides Viral Clearance and TSE Clearance studies, Mycoplasma, DNA and Bacterial Clearance studies can also be performed.




Viral Clearance Study Design



Viral Clearance Studies

When are they performed?



Viral Clearance Studies

Testing each process step

Spike virus into process intermediate. Must be appropriate to Load material (e.g., liquid vs. solid)

Input Virus

Cleared Virus

Leftover Virus

Perform scaled-down chromatography process step(s)



Collect fractions (e.g., flow-through, eluate, regenerate)



Solvent, Detergent, pH inactivation




Perform scaled-down filtration process step

Collect fractions

Titrate for virus infectivity level (TCID50)*



$$\text{Reduction} = \frac{[\text{virus}] \text{ load}}{[\text{virus}] \text{ product}}$$



Viral Clearance Studies

Generic timelines


Study Type	Wk 1	Wk 2	Wk 3	Wk 4	Wk 5	Wk 6	Wk 7	Wk 8	Wk 9	Wk 9	Wk 10	Wk 11	Wk 12	Wk 13	Wk 14	Wk 15	Wk 16	Wk 17	Wk 18	Wk 19	Wk 20
------------	------	------	------	------	------	------	------	------	------	------	-------	-------	-------	-------	-------	-------	-------	-------	-------	-------	-------



Conclusion

- Three-Pillar Approach for Total Virus Safety Considerations in Plasma
- Start with the end in mind
- Merck Key Competencies in Plasma:

Regulatory Experts, Virus Clearance Team, Process Validation Experience, Plasma Fractionation Process Knowledge and Leading Technologies





Thank YOU!

contact:
Louis.Wong@Merckgroup.com



MERCK